

Axial (+3) Image



Axial (+3) Ground Truth



Axial (+3) Eval Composite



Patient: CQ500_0030
Dice: 0.6026 | IoU: 0.4312
Prec: 0.9001 | Rec: 0.4529 | Spec: 0.9999
TP:2640 FP:293 FN:3189 TN:4188182

GT Vol: 15.03 mL | Pred Vol: 7.56 mL
RVD: -49.7%
HD95: 91.03 mm | ASSD: 10.57 mm

Lesions GT/Pred: 1/2
Lesion Sens/Prec: 1.00/0.50

Optimal Thresh: 0.30 (Dice 0.6400)

Diagnosis:
Severe under-segmentation; missed significant hemorrhage volume.
Boundary distance is large (HD95=91.0mm), indicating localized regions of substantial under-segment even where bulk overlap looks fine.

- True Positive
- False Positive
- False Negative
- Ground Truth

Axial (Center) Image



Axial (Center) Ground Truth



Axial (Center) Eval Composite



Axial (-3) Image



Axial (-3) Ground Truth



Axial (-3) Eval Composite

